


Photongreen

EPILAZIONE NATURALE

Technical Data

- **Type:** LED emitter
- **LED wavelength:** 525 Nm
- **Pulsation:** approx. 1 sec.
- **Dimensions:** 27x32x100 cm
- **Weight:** 27 kg
- **Power:** 220 Volt A.C.
- **Absorbed power:** 30 Watt max



Why should I have Photongreen in my beauty centre?

Photongreen is the equipment that prevents and delays the regrowth of hair, without causing any damage.

RECOMMENDED FOR:

- Any skin type and colour
- Any hair type and colour
- Helps against reddening
- Helps against couperose


What technology is used?

Photongreen technology is based on the light that prevents the regrowth of hair.


Due to this revolutionary treatment, the hair no longer receives the nutrition it needs to grow and it therefore weakens to the point at which it doesn't grow back for a long time.

LIGHT → SCLEROTIZATION HAIR

Nutrition of the follicle through the blood vessels



Interrupting the nutrition of the follicle



L = Lymphatic - A = Arterial - V = Venous. Before and after.

How is it used?

After waxing*, due to the joint action of the photo-carrier lotion (PhotonBlue) and the pulsating flashes of green light at 525 Nm, Photongreen delays the regrowth of hair, without causing any damage or burns, thanks to the hair sclerotization process

*The Photongreen action depends on the waxing having been done perfectly. In fact there must be absolutely no hair in the follicle at the time of application.

Wax perfectly



Use Photongreen




How long is needed for each treatment?

Average treatment time recommended per zone:

- Leg = 1/2 minutes approx.
- Groin = 30 seconds approx.
- Armpit = 30 seconds approx.

STAGES OF HAIR



PhotonBlue lotion

The properties of the photo-carrier lotion make it an essential tool for the application of green light to prevent and delay the regrowth of hair.

The lotion actually carries the light and the substances contained in it inside the follicle.

Here are the main properties:

- Good luminous conductivity
- Dermatologically tested
- Contents of substances of vegetable origins
- Hypoallergenic
- Perfume-free

Photongreen also offers a Blue light mouse for dark skin or facial hair.

Testimonies

I'm Denise Bertolazzi, from the **DENISE** Beauty Centre in Timoline Corte Franca (BS).

I've had excellent results with Photongreen over just a few months: thinner hair, slow regrowth, to my complete satisfaction and that of my customers: a real success! I've been able to quickly repay the investment and it's all earnings from now on.

THANK YOU, PHOTONLIFE!

Denise Bertolazzi

Denise Bertolazzi
DENISE Beauty Centre, Timoline Corte Franca (BS)

PhotonLife promotes you

PhotonLife works! So that you can tell your end customers all about it, it's essential that you have the right graphic-based publicity in the customer's beauty centre. Here is a sample of the **basic bespoke kit** to effectively display the results that are achieved.



Explanatory flyers
F.to 20x20 cm



Roll-up
F.to 80x200 cm



Counter stand
F.to 21x30 cm



Display case stand
F.to 60x90 cm

PhotonLife at your service

There's no such thing as **success without effective communication**, especially in a sector like beauty. Clear, professional and targeted initiatives create the fundamental basis to engage end customers in beauty centres, causing them to test the benefits of the services on offer. This is why, as well as the state-of-the-art equipment, **we will provide assistance in:**




Professional training

To use the devices correctly and for a better in-depth knowledge of their technology.



Open day demonstrations

Days, events with demonstrations in your centre to win over customers, old or new.



After-sale service

Telephone marketing assistance and free technical advice always at your disposal.



Directed advertising campaigns

Editorials in the main women's and men's weeklies, creating a demand, national advertising campaigns.



PhotonLife Club

Only a "team" of trained and capable professionals can feel up to explaining the exceptional possibilities of the PhotonLife equipment. The idea of a **CLUB** came about, bringing you the following advantages:

- **Certifications**
- **Courses**
- **Events**
- **Special discounted conferences**
- **Directed advertising campaigns**
- **The chance to use the PhotonLife brand in your beauty centre**

www.photonlife.it

PhotonLife®

LABLIFE srl - Via Pieve, 5 - 25080 Nuvolento (BS) - Italy - Tel. +39. 030.6898 666 - Fax +39. 030. 6898 054 - info@lablife.it



Photongreen
EPILAZIONE NATURALE

SALES PROMOTION